


LINSEIS

T H E R M A L A N A L Y S I S

NUCLEAR

DIL
STA
LFA



Since 1957 LINSEIS Corporation has been delivering outstanding service, know how and leading innovative products in the field of thermal analysis and thermo physical properties.

We are driven by innovation and customer satisfaction.

Customer satisfaction, innovation, flexibility and high quality are what LINSEIS represents. Thanks to these fundamentals our company enjoys an exceptional reputation among the leading scientific and industrial organizations. LINSEIS has been offering highly innovative benchmark products for many years.

The LINSEIS business unit of thermal analysis is involved in the complete range of thermo analytical equipment for R&D as well as quality control. We support applications in sectors such as polymers, chemical industry, inorganic building materials and environmental analytics. In addition, thermo physical properties of solids, liquids and melts can be analyzed.

LINSEIS provides technological leadership. We develop and manufacture thermo analytic and thermo physical testing equipment to the highest standards and precision. Due to our innovative drive and precision, we are a leading manufacturer of thermal Analysis equipment.

The development of thermo analytical testing machines requires significant research and a high degree of precision. LINSEIS Corp. invests in this research to the benefit of our customers.




Claus Linseis
Managing Director



German engineering

The strive for the best due diligence and accountability is part of our DNA. Our history is affected by German engineering and strict quality control.



Innovation

We want to deliver the latest and best technology for our customers. LINSEIS continues to innovate and enhance our existing thermal analyzers. Our goal is constantly develop new technologies to enable continued discovery in Science.



THERMAL ANALYSIS OF NUCLEAR MATERIALS

Special field of application: Nuclear materials

Since the 1950s nuclear energy is the world-wide most important energy source around the world. With their advantage of clean and cheap power supply, core reactors were undergoing a continuous global improvement during the last 50 years. Meanwhile the reactors of the 4th generation such as very high temperature reactors (VHTR) or sodium cooled fast reactors (SFR) as well as the unique molten salt reactor (MSR) are currently under development and will be the fu-

ture for nuclear energy.

Due to the research that is done in that field, there is a need of analytical equipment and especially instruments for thermal analysis. Of course these special applications and safety requirements need a lot of modifications of the standard devices, that makes Linseis become the worldwide leader in thermal analysis of nuclear materials as we are the most flexible and most experienced player on that market.

LFA




Standard Test Method for Thermal Diffusivity by the Flash Method

A small, thin disc specimen is subjected to a high intensity short duration radiant energy pulse. The energy of the pulse is absorbed on the front surface of the specimen and the resulting rear face temperature rise is recorded. The thermal diffusivity value is calculated from the specimen thickness and the time required for

the rear face temperature rise to reach certain percentages of its maximum value.

Out of density and Cp the thermal conductivity can be calculated using the following formula:

$$\lambda(T) = a(T) \cdot c_p(T) \cdot \rho(T)$$




Calculation of thermal diffusivity

- Determine the baseline and maximum rise to give the temperature difference, ΔT_{\max}
- Determine the time required from the initiation of the pulse for the rear face temperature to reach $\Delta T_{1/2}$. This is the half time, $t_{1/2}$.
- Calculate the thermal diffusivity, a , from the specimen thickness, L squared and the half time $t_{1/2}$, as follows:

$$a = 0.13879 L^2/t_{1/2}$$


Combined solution of the simultaneous heat loss and finite pulse corrections with the laser flash method



Conclusion

The combined model method with nonlinear parameter estimation has been proven for more than 100 samples. In all cases it worked reliably and its results gave the correct adiabatic, finite pulse, and/or heat loss corrected values. The


two main advantages of the method are that no operator choice between the different models and correction is necessary, and the fit can be checked by plotting the model curve.



Calculation of thermal diffusivity

- Determine the baseline and maximum rise to give the temperature difference, ΔT_{\max}
- Determine the time required from the initiation of the pulse for the rear face temperature to reach $\Delta T_{1/2}$. This is the half time, $t_{1/2}$.
- Calculate the thermal diffusivity, a , from the specimen thickness, L squared and the half time $t_{1/2}$, as follows:

$$\alpha = 0.13879 L^2/t_{1/2}$$



For Special setup for radioactive or toxic samples


- Measurement unit separated from control electronics and laser glovebox
- Laser connected by fibre glass optical cable
- Measurement unit can be placed in hood or glovebox
- Maintenance, operation and setup possible in glovebox

TECHNICAL SPECIFICATIONS

	LFA 500 IR	LFA 1000	LFA 1000 HT
Sample dimensions	<p>∅ 10 / 12.7 / 25.4 mm, 0.1 to 6 mm thick</p> <p>10 / 20 mm; 0.1 up to 6mm thick</p>		
Samples	solids, liquids, powders. pastes	solids, liquids, powders. pastes	solids, liquids, powders. pastes
Sample robot	up to 6 samples	up to 6 samples	up to 3 samples
Vacuum	depends on model up to 10 ⁻⁴ mbar		
Atmosphere	inert, oxidizing or reducing		
Measuring range Thermal Diffusivity	0.01 up to 1000 mm ² /s	0.01 up to 1000 mm ² /s	0.01 up to 1000 mm ² /s
Measuring range Thermal Conductivity	0.1 to 2000 W/(m·K)	0.1 to 2000 W/(m·K)	0.1 to 2000 W/(m·K)
Pulse source	Flash lamp	Flash lamp	Flash lamp
Pulse energy	10 J/pulse	25 J/pulse	25 J/pulse
Pulse energy & pulse duration adjustment	yes	yes	yes
Pulse length adjustment	software adjustable	software adjustable	software adjustable
Sensor type	InSb / MCT	InSb / MCT	InSb / MCT
Furnace model	<p>IR-furnace RT up to 500°C RT up to 1000°C</p> <p>Cryo-Furnace -125 up to 600°C</p>	<p>Resistance heater RT up to 1250°C RT up to 1600°C</p>	<p>Graphite RT up to 2000°C RT up to 2800°C</p>

APPLICATIONS


Graphite standard



The curve shows the thermal diffusivity data of a graphite standard from NIST that was measured with a special LFA with fibre-optical connected laser and separated electronics. The results are matching with the literature values and accuracy and power of the laser are completely identical with the standard instrument.

Graphite is one of the most important materials in reactor construction and due to its high thermal conductivity and temperature stability it offers various purposes it can be used for.

Thermal Conductivity of three carbide ceramics



Three different types of carbides were measured by LFA to determine the thermal conductivity. One of the samples has been the uranium derivative that showed a considerable low value of around 25 W/mK. ZrC shows an increasing trend that has an electronic dominant character, whereas the SiC shows a more or less common decrease shape.

DILATOMETRY



Technical Expansion Coefficient

The technical expansion coefficient can be calculated as follows:

$$\alpha_{\text{tech}}(k) = \frac{1}{L_0} \frac{\Delta L_k - \Delta L_0}{T_k - T_0} \quad (k=1, \dots, n)$$

for every measured point k.

L_0 is the sample length at 20°C, ΔL_0 the change in length at 20°C (linear extrapolated out of the first data points), ΔL_k the according length change at the temperature T_k .

Measuring thermal expansion of radioactive samples

The standard instrument for measuring the thermal expansion coefficient (CTE), the Dilatometer, detects the change in length of any sample by LVDT detector systems. For toxic and radioactive samples, it was necessary to remove all electronics and build a special furnace that is separated from the detector. Beside this, the design was also changed to a more accessible

setup for maintenance and mechanical changes under glovebox or fume hood.

Another method for CTE determination is the contact free optical dilatometer. It can be used with any atmosphere and the sample does not see any force. For less x-ray emitting samples it can be a very useful instrument, however it can also be placed within a glovebox or fume hood.

Measuring systems

The LINSEIS L75 Dilatometer can be used with different measurement systems for various applications, even for measuring powders and pastes, offering a broad range in temperature and field of application.



Quartz measuring system

ø 7/12 mm



Al₂O₃ measuring system contact

free



Al₂O₃ measuring system standard



Adapter for powders and pastes



Quartz measuring system for


large samples ø 20mm

TECHNICAL SPECIFICATIONS

	L 75 Horizontal	L75 Vertical	L75 Laser	L74 Optical
Temperature range	-180 up to 2000°C	-180 up to 2800°C	-180 up to 1000°C	-100 up to 2000°C
Atmospheres	Inert, oxidizing, reducing	Inert, oxidizing, reducing, vacuum	Inert, oxidizing, reducing, vacuum	Inert, oxidizing, reducing, vacuum
Resolution	1nm	1nm	0.3nm	1µm

Thermal expansion Measurements

We find that the Coefficient of Thermal Expansion (CTE) of super-invar is sensitive to the level of irradiation exposure.




Customer Application from Brookhaven:


The CTE of invar is dependent on the level of irradiation exposure.

The experiments and studies were performed using our L75 horizontal as a special version in a glovebox.

APPLICATIONS

Super-invar







STA



DSC- True Heat Flow measurement

Quantitative DSC-signal



Differential Scanning Calorimetry (DSC)

"A technique in which the difference in energy input into a substance and a reference material is measured as a function of temperature, while the substance and reference material are subjected to a controlled temperature program."

Differential Signal

The differential signal is displayed as a baseline. Effects, for example the melting of a metal, can be observed as a peak. The area of the peak gives the amount of enthalpy and the direction of the peak indicates the way of heat flux – endothermic (down) or exothermic (up).

Temperature vs. Time

During an effect like a reaction, decomposition or phase transition, a temperature difference (heat flux difference) between the sample and the reference crucible can be measured by means of a thermocouple.

MEASURABLE PROPERTIES

- Mass change as % and mg
- Rate controlled mass loss
- Evaluation of mass loss
- Residue mass evaluation
- Compositional analysis



- Enthalpy
- Endo- / Exo- thermic
- Phase transformation
- Melting point
- Glass point
- Crystallinity


- Thermal stability
- Oxidation stability
- Purity
- Solidus / Liquidus relationship
- Product identification

SENSORS


Our STA can be equipped with an unmatched amount of different user exchangeable TG-DSC, TG-DTA or TG sensors.


Each sensor is available with different thermocouples to provide the highest sensitivity for your desired temperature range.





Best possible sensitivity for your application






TECHNICAL SPECIFICATIONS


	STA PT 1000	STA PT 1600			STA HP/2	
Temperature range	RT up to 1000°C	-150 up to 2400°C			RT up to 1100/ 1400°C/1600°C/1800°C	
Vacuum	optional 10 ⁻² mbar	10 ⁻⁵ mbar (depends on vacuum pump)			up to 10 ⁻⁴ mbar	
Pressure		up to 5 bar (optional)			up to 150 bar custom solution on request	
Heating rate	0.01 up to 100°C/min	0.01 up to 100°C/min (depends on furnace)			0.01 up to 100°C/min (depends on furnace)	
Temperature precision	0.01°C	0.01°C			0.05°C	
Sample robot	optional 42 / 84	optional 42			-	
TG						
Resolution	0.1 µg	0.025 µg	0.1 µg	0.1 µg	0.1 µg	0.1 µg
Sample weight	5 g	5 g	25 g	35 / 50 g	15 / 25 g	35 / 50 g
Measuring range	25 / 2500 mg	25 / 2500 mg	25 / 2500 mg	35000 mg	25 / 2500 mg	35000 mg
DSC						
DSC-sensors	E / K / S	E / K / S / B / C			E / K / S / B / C	
DSC resolution	0.3 / 0.4 / 1µW	0.3 / 0.4 / 1 / 1.2 µW			0.3 / 0.4 / 1 / 1.2 µW	
Calorimetric sensitivity	approx. 4 / 6 / 17.6 µW	approx. 4 / 6 / 17.6 / 22.5 µW			approx. 4 / 6 / 17.6 / 22.5 µW	
DTA						
DTA-resolution	0.05 µV	0.05 µV			0.05 µV	
Sensitivity	1.5 µV/mW	1.5 µV/mW			1.5 µV/mW	
DTA-measuring ranges	250 / 2500 µV	250 / 2500 µV			250 / 2500 µV	

INCONEL 718




The curve shows the specific heat of Inconel 718 during heating. The steel can be used to produce super invar steel that can be hardened by irradiation. Out of the specific heat curve, the solidus and liquidus temperature at 1254°C and 1343°C as well as the melting enthalpy of 186J/g can be determined.

Cement



The main parts of cement are tri calcium silicate, di calcium silicate and tri calcium aluminates. After putting on the cement with water different hydrates slowly form. The absorbed water evaporates first, then hydrates of the calcium silicate decompose and at 570°C the hydroxides of calcium, magnesium and aluminum follow. Subsequently, calcium carbonate CO₂ splits off.

Decomposition of gypsum and quartz sand



Typical Applications for Gypsum and quartz sand are plaster and mortar. The content of gypsum in the sample illustrates a two-step release of H₂O from CaSO₄ · 2H₂O (dihydrate) into CaSO₄ · ½H₂O (halfhydrate) and finally into CaSO₄ (anhydrite). This requires an entire energy of 121.6 J/g. The Quantitative analysis (TG-DSC Sensor) reveals that the sample contained 23.4% of pure dihydrate. Between approx. 300°C and 450°C, the exothermic formation of β-CaSO₄ with a released energy of 18.32 J/g occurred. The endothermic effect at a temperature of 574.9°C is due to the structural α → β transition of quartz (crystalline SiO₂).

**LINSEIS GmbH Germany**

Vielitzerstr. 43

95100 Selb

Tel.: (+49) 9287 880 0

E-mail: info@linseis.de

**LINSEIS Inc. USA**

109 North Gold Drive

Robbinsville, NJ 08691

Tel.: (+1) 609 223 2070

E-mail: info@linseis.de

**LINSEIS China**

Kaige Scientific Park 2653 Hunan Road

201315 Shanghai

Tel.: (+86) 61 90 12 03

Tel.: (+86) 50 55 06 42

E-mail: info@linseis.com.cn

**LINSEIS France**

1 Route de Trévoux

69250 Neuville/Saone

Tel.: (+33) 6.24.72.33.31

E-mail: contact@ribori-instrumentation.com

**LINSEIS Poland**

ul. Dabrowskiego 1

05-800 Pruszków

Tel.: (+48) 692 773 795

E-mail: info@linseis.de

www.linseis.com**Products:** DIL, TG, STA, DSC, HDSC, DTA, TMA, MS/FTIR, In-Situ EGA, Laser Flash, Seebeck Effect, Thin Film Analyzer, Hall-Effect**Services:** Service Lab, Calibration Service

08/17

LINSEIS